

---

# Software Release Notes

---

## CEMS Portal Version 5.0.0

This document provides the changes to CEMS version 5.0.0.

### Purpose of the Release

The purpose of this document is to communicate major new features, changes and solved bugs in this release of CEMS software. It may also document workarounds and any feature restrictions.

### Software Updates

- Change originator from serial number to device name
- Update SMS messages to fit 160-character limit
- Rename of roles used to identify different administrators and users of site
- Various security issues were identified and repaired
- Rename of titles on tiles or menus
- Dynamic drag and drop list on floor plan

### New Software Features

- Calibration Tracking Tool
  - Ability to track, alarm and report calibration on and through dates for devices
- Enable Integrator Partnerships
  - Ability to easily setup and manage multiple customer dashboards
- Addition of new menu tabs
  - General, Users, Calibration Tracking Tool, Help
- Addition of new items on Floor Plan legend
- Ability to add user defined devices

### Solved Software Bugs

- Issue related to Unit of measure as “Network” and alarm status in reporting
- For Flex Device, CEMS doesn’t match EDGE in terms of user defined parameters
- Issues related to halo dot
- Various minor issues were identified and repaired

Setra Systems, Inc.

159 Swanson Road, Boxborough, MA 01719  
800.257.3872 • www.setra.com