
Software Release Notes

CEMS Portal Version 5.0.2

This document provides the changes to CEMS version 5.0.2.

Purpose of the Release

The purpose of this document is to communicate new features, changes, and solved bugs in this release of CEMS software. It may also document workarounds and any feature restrictions.

Software Updates

- Default Check boxes for alarms table filter

New Software Features

- None

Solved Software Bugs

- Report Scheduler: XLSX reports are inconsistent
- Unreadable text gets added to SMS text messages using AT&T gateway for alarm/action notifications
- Inactivity timeout is not correctly mapped with the time gap in individual trend charts
- Data exported from a parameter is missing data and /or incorrectly ordered
- In a session, for PDF/DOCX reports for user defined device types, the maximum and average values are missing in the data summary for User Defined 7 parameter
- If 0 is entered for an alarm and/or action level threshold, the line in the trend chart vanishes
- Report scheduler triggered CSV md5sum fails while validating through report validator