

# SOFTWARE RELEASE NOTES

## CEMS Portal Version 6.0.1

This document provides the changes to CEMS version 6.0.1

### Purpose of the Release

The purpose of this document is to communicate solved bugs in this release of CEMS software. It may also document workarounds and any feature restrictions.

### Solved Bugs

- Mistype for Reasons Pre-Configuration Menu
- Pop-up notifications from devices can be viewed across all campuses
- Audit log sign-off is not properly recording the reasons for alarm events
- Clean up of reminiscence of action level to warning naming convention
- Non-clearable fault alarm indication in tile view and halo dots

Setra Systems, Inc.

159 Swanson Road, Boxborough, MA 01719

800.257.3872 • [www.setra.com](http://www.setra.com)