



MRMS BACnet Object List

Object Type	Instances	Object Name	Object Description
Device	DEV1	Device	Description of the Nurse's Unit
Analog Value (34 bytes/instance)	AV1	Audible Alarm	The buzzer volume. It can be changed from 0 to 4, where '0' disables the Alarm sound and values 1, 2, 3 & 4 indicate the buzzer volume of the nursing station device.
	AV2	Mute Time Out	Buzzer volume mute timeout (Ranges from 0 to 9999). 9999 will make the buzzer volume mute forever.
	AV3	Number of Rooms	The number of rooms configured. Maximum can be up to 8 rooms.
Binary Value (16 bytes/instance)	BV1	Supervisor Password	0 – Disabled 1 - Enabled
	BV2	Operator Password	0 – Disabled 1 - Enabled
Character String Value/Proprietary Object	PCV1	Unit Label	
Assuming (100 bytes/instance)	PCV2	Technical Support	