



### AOLM48 Series - Olympus Microscopes

Mounting: Threaded Locking Ring (Included) Working Distance: 1.25in. – 5.0in.

Model	Part #	Cable Length (in.)	End Tip Dia. (in.)	Light Source Adapter .590 (15mm) Input	Light Source Adapter 1.00 (25mm) Input	Applications
AOLM48	700008440048	48	.588	Direct Fit	SX-13	SZ-30, SZ-40, SZ-60

### LP25Li36M Series - Leica Microscopes

Mounting: 1 Thumb Screw Working Distance: 1.25in. - 5.0in.


Model	Part #	Cable Length (in.)	End Tip Dia. (in.)	Light Source Adapter .590 (15mm) Input	Light Source Adapter 1.00 (25mm) Input	Applications
LP25Li36M	700014060036	36	.588	Direct Fit	SX-13	S4E, S6D/E, S8APO, GZ-6 (.5X and 2.0X)
LP25Li72M	700014060072	72	.588	Direct Fit	SX-13	S4E, S6D/E, S8APO, GZ-6 (.5X and 2.0X)

### AUNi48 Series - Unitron Microscopes


Mounting: Threaded Locking Ring (Included) Working Distance: 1.25in. – 5.0in.

Model	Part #	Cable Length (in.)	End Tip Dia. (in.)	Light Source Adapter .590 (15mm) Input	Light Source Adapter 1.00 (25mm) Input	Applications
AUNi48	700008442048	48	.588	Direct Fit	SX-13	ZST, ZSB


# AOLM48 Dimensions:



# LP25Li36M Dimensions:



# AUNi48 Dimensions:



**Note:** Ring lights available in variety of lengths. 36 inch version shown in diagram above. Please contact Dolan-Jenner for availability on length variations. Please note two standard end ferrules available along with many custom designs.

**Dolan-Jenner**  
industries

Dolan-Jenner Industries  
159 Swanson Rd.  
Boxborough, MA 01719 USA

Phone: 800.833.4237  
Fax: 978.264.0292  
Email: support@dolan-jenner.com

1.800.83.FIBER

WWW.DOLAN-JENNER.COM